

Product safety Consumer satisfaction Pharmaceutical purity Accurate auditing and connectivity Sanitary processing Strain typifying

Laboratory and Pharmacy





To satisfy requirements for selling cannabis

derivatives, edibles and topicals by weight, Rice Lake laboratory balances can weigh precisely to 0.001 gram. Data connectivity for labeling and tracking compliance is available with RS-232, Bluetooth® or USB port.

Retail Shops and Dispensaries

Balances Point of Sale Scales

-il ----l--



Rice Lake retail scales are ideal for the point of

sale (POS) or retail packaging counters in cannabis storefronts and dispensaries. Legal for Trade scales allow legal selling of products by weight from the scale, for the highest level of accuracy in weight-based transactions. Rice Lake's retail scales have data communication options suitable for tracking, reporting and inventory management, and select models o er additional weight displays for transparency and customer con dence.





Verification and Compliance



Castiron Calibration Weights Precision Calibration Weights Leaf Calibration Weights

Rice Loke offers a full line of calibration weights and services required for keeping scales compliant as well as for best practices in sales and quality assurance. With weights available in almost any class, set, size or individual weight, Rice Lake's knowledgeable RYERY Graubestime weight of th

scale verification need.

Unique Production Needs RESOLVED

For thefastest returnoninvestmentincannabisproduction facilities, turn toequipmentthattrulyaids profitability—accurate and durable weighing equipment designedfor yourfacility'senvironment,fullybacked by Rice Lake's legendary quality and accessible, comprehensive customersupport.

Explorer® Series, OHAUS Precision Balance





Standard Features

- · Large, 5.7 full-color touchscreen display with intuitive setup and application software
- · AutoCal[®] automatic internal adjustment weight for optimized weighing performance
- 14 application modes, 11 weighing units and 5 display languages for a wide variety of applications
- · Four programmable touchless sensors for hands-free operation
- 0.001 g models feature frameless draftshield with anti-static glass, flip-top and large side entry doors on top-mounted bearings

Options

· External calibration models

FX-i Series, A&D Weighing Precision Balance



As of: Wednesday, May 28, 2025



Standard Features

- Compact Super Hybrid Sensor Technology (SHS)
- Vacuum Fluorescent Display (VFD)
- GLP/GMP/ISO compliant to either AD-8121 printer or computer. Transmits date, time, ID number, serial number and calibration data.

Approvals



RS-232 interface for connecting to printer or PC

- Response indicator function
- Weighing chamber
- · Underweigh hook

Options

- Built-in rechargeable battery
- Quick USB interface (FXi-02) option
- LAN-Ethernet interface
- WinCT© communication tools software
- Note: FXi NTEP and FXi models available by special order (internal adjustment weight)





JOSE RAFAEL RAMIREZ

GENERAL MANAGER

§ 516 927 4705
⊕ www.vivatechrs.com/

516 413 7006jose.ramirez@vivatechrs.com